




Application Note

Introduction

The aim of our experiment is to evaluate antibody dependent phagocytosis of SARS-CoV-2 by monocytes. As a first step, we had to check whether the IgA PTXCOV-A569 binds to the S1 spike and at which concentration in ELISA.

Results

As the graph below displays, the IgA PTXCOV-A569 showed a good signal at $1\mu g/ml$ at OD_{450nm} after 15 min of reaction development.



IgA-anti-S1: Anti-SARS-CoV-2 S1 antibody, Human IgA (PTXCOV-A569) using at 1ug/ml and 10ug/ml Anti-SARS-CoV-2 ELISA IgA (Euroimmun kit El2606-9601A) following manufacturer's instructions

Method

We tested the IgA PTXCOV-A569 in ELISA using a kit. The plate was precoated with SARS-Cov-2 S1 spike protein. We added the IgA PTXCOV-A569 at 1μ g/ml and 10μ l/ml.

We used the reagents from the kit to detect binding of IgA PTXCOV-A569. We used the positive and negative controls from the kit as well. Binding was evaluated by spectrophotometry at OD_{450nm}.

Conclusion

We used the entirety of the free sample provided by **ProteoGenix** for this ELISA, we will buy more of the product in order to pursue our experiments and reach our initial aim.

Once we have more IgA PTXCOV-A569, we will continue to test the IgA PTXCOV-A569 at 1μ g/ml in our function assay of antibody dependent phagocytosis.



